



# CAREONE DIAGNOSTICS

PRICE-LIST w.e.f 01/05/2013

The collection of all tests and their prices and methods of collection of samples.

**Careone**  
**1/5/2013**

# CAREONE DIAGNOSTICS

## SPECIAL TEST PRICE-LIST

Test Code No	TEST NAME	RATES IN RUPEES	SPECIMEN	Specimen Vol	TEMP	SAMPLE INSTRUCTION	REPORTING TIME
1	1,25 DIHYDROXY CHOLECALCIFEROL (VITAMINE D)	3500	Serum	2ml	Rf / Fro		10 Days
2	17- HYDROXYPROGESTRONE; 17- OHP Serum	900	Serum	2ml	Rf / Fro		4 days
3	17-HYDROXYCORTICOSTEROIDS , 24 hrs urine	6000	Urine	50ml	Ambient	50 ml (20 ml min.) aliquot ,24 hrs urine collected in 50% HCL	7 days
4	17-KETOSTEROIDS, 24 HRS URINE	1600	Urine	50ml	Ambient	50 ml (20 ml min.) aliquot ,24 hrs urine collected in 50% HCL	One week
5	5- ALPHA DHT, DI HYDRO TESTOSTERONE, SERUM	1200	Serum	2ml	Rf / Fro		One week
11	ACE, ANGIOTESIN CONVERTING ENZYME	700	Serum	2ml	Rf/Fro		Next Day
12	ACETAMINOPHEN (PARACETAMOL)	1000	WB-EDTA	3ml	Rf/Fro		One week
14	ACETYL CHOLINE RECEPTOR ANTOBODIES	5000	Serum	2ml	Rf/Fro		One week
18	ACTH Stimulation Test for Cortisol	350 (each)	Serum	2ml	Ambient	Baseline sample Administer (Synacthen) 250 ug IV,draw additional samples at 30 and 60 min.	Same day
19	ACTH, PLASMA	1000	Plasma-EDTA	2ml	Frozen	Draw blood between 8-10 a.m. into 1 Lavender Top (EDTA) tube. Mix well, separate and transfer plasma IMMEDIATELY into a plastic transport container and FREEZE. Ship frozen. DO NOT THAW.	3 days
21	ACTIVATED PROTEIN C-RESISTANCE; APC RESISTANCE ASSAY; APCR Factor V leiden (FD), FACTOR V LEIDEN (Functional assay)	3500	Citrated - Plasma in 2vials	4ml	Frozen	Platelet Poor Plasma (PPP) from 2 Blue Top (Sodium Citrate) tubes prepared within 1 hour of collection.	7 days
22	ADH- ANTIDIURETIC HORMONE (VASOPRESSIN)	4000	Plasma-EDTA	3ml	Frozen		8 days
23	ADENOSINE DEAMINASE ; ADA, SERUM/FLUID	400	WB-EDTA or Serum,Fluid	2ml/2ml/5ml	Frozen	Hemolysed sample not accepted, 5 ml fluid in a sterile container.	Same day
24	ADRENALINE (EPINEPHRINE) PART OF	3000	WB-EDTA	3ML	Frozen	Separate plasma freeze and transport	5 days

# CAREONE DIAGNOSTICS

## SPECIAL TEST PRICE-LIST

	CATECHOLAMINE					frozen.	
29	ALCOHOL SCREEN, BLOOD	500	WB-SodiumFluride	3ml	Ambient	Cleanse collection site with NON-ALCOHOLIC agent before drawing blood. Ship refrigerated. DO NOT FREEZE.	Next day
30	ALCOHOL SCREEN, RANDOM URINE	500	Urine	10ml	Rf	Aliquot of a random urine sample filled up to the brim of a sterile screw capped container. DO NOT FREEZE.	Next Day
31	ALDOSTERONE, SERUM	1300	Serum	2ml	Rf/Fro	2 mL ( 0.5 mL min.) serum from 1 SST. Overnight fasting is preferred. Draw upright specimen after the patient is made to sit up for 30 min. Ship refrigerated or frozen.	One week
32	ALDOSTERONE, 24 HRS, URINE	1500	Urine	5ml	Rf/Fro	5 ml aliquot of 24 hrs urine refrigeratde during collection. Add 1 gm boric acid per 100 ml of urine. Mix well, pH 7.5 mix well and shed refrigerated or frozen. Mention 24 hrs volume on test request form and on the container	7 days
34	ALKAPTONURIA/ HOMOGENITISIC ACID, URINE	300	Urine	10ml	Ambient	Aliquot of freshly voided urine in a plastic container . No preservative, Clinical and drug history mandatory	2 days
35	Allergy individual markers Acacia, Aspergillus, Bermuda Grass Cat Epithelium, candida, chicken, Cladosporium Cockroach, Cow dander Dog Dander, Egg white Eucalyptus Horse dander, Sorgum, milk, mites, Pea nut, Rice Shrimp, soyabean, wheat	900 each	Serum	3ml	Rf/Fro		One week
36	Allergy Screen Adult	6000	Serum	3ml	Rf/Fro		One week
37	Allergy Screen paediatrics	6000	Serum	3ml	Rf/Fro		One week
38	Food Panel	2500	Serum	3ml	Rf/Fro		One week
39	Animal panel	2500	Serum	3ml	Rf/Fro		One week
40	Mould Panel	3000	Serum	3ml	Rf/Fro		One week

# CAREONE DIAGNOSTICS

## SPECIAL TEST PRICE-LIST

41	Grass Panel	2500	Serum	3ml	Rf/Fro		One week
42	Tree Panel	2500	Serum	3ml	Rf/Fro		One week
43	ALPHA - FETOPROTEIN, AMNIOTIC FLUID; AFP, AMNIOTIC FLUID	500	Amiotic Fluid	2ml	Rf/Fro	Avoid contamina- tion with maternal blood.	Same day
44	ALPHA - FETOPROTEIN, SERUM AFP Cancer Marker	500	Serum	3ml	Rf/Fro	Give brief clinical history.	Same day
45	ALPHA - FETOPROTEIN, SERUM AFP Maternal Screening.	500	Serum	3ml	Rf/Fro	Give brief clinical history.	Same day
46	ALPHA-1-ANTITRYPSIN QUANTITATION; AAT	1500	Serum	3ml	Rf/Fro	2 mL (0.5 mL min.) serum from 1 Red Top (No Additive) tube. Overnight fasting is preferred. Spin and separate serum immediately. Ship refrigerated or frozen.	Next Day
47	Aluminium, dialysis fluid	1500	Serum	10ml	Rf/Fro	10 ml (5 ml min) dialysis fluid in a special royal blue top tube available from lab	10 Days
48	Aluminium, plasma/ SERUM	2000	Serum/Hep arin	3ml	Rf/Fro	In special Royal blue top ( HEPARIN) tube for plasma or 1 royal blue top ( NO ADDITIVE) tube for serum available from lab	10 Days
49	AMMONIA, BLOOD	500	WB-EDTA	3ml	Refrigerat ed	Avoid haemolysis. KEEP ON ICE. DO NOT FREEZE.	Same day
50	AMMONIA	500	Plasma- EDTA	3ml	Ambient		Same day
52	Amphetamine and metabolite screen, random urine	500	Urine	10ml	Rf/Fro	10 ml aliquot of random urine.	Next Day
55	AMYLOID STAIN (CONGO RED)	500	Tissue,4Gla ss Slides,Parra fin Block	4	Ambient	Submit tissue biopsy in 10% formal saline OR formalin fixed paraffin embedded tissue block OR FNAC specimen smeared onto 4 glass slides and air dried. Ship at room temperature. Brief clinical history in Histopathology Requisition Form is mandatory.	3 Days
56	ANCA BY IMMUNOFLUROSCENCE ( anti neutrophillic cytoplasmic antibody) (Incl. P- ANCA & C-ANCA)	1500	Serum	4ml	Rf/Fro	4 ml gel ( 2 ml min. ) serum	Sample by Tue/Friday Reporting Wed/Sat
57	ANDROSTENEDIONE	1000	Serum	3ml	Rf/Fro	An early morning specimen is preferred.Specify age and sex on test request form.	3 days

# CAREONE DIAGNOSTICS

## SPECIAL TEST PRICE-LIST

59	ANTI-BODIES TO EXTRACTABLE NUCLEAR ANTIGENS, COMPREHENSIVE PANEL; ENA COMPREHENSIVE PANEL Individual antibodies to * Sm *RNP-Sm Complex * SS-A/Ro * SS-B/La *Jo-1 * Scl-70 * Centromere	5000	Serum	2ml	Rf/Fro	Overnight fasting is preferred.	7 Days
60	ANTI-BODY CYCLIC CITRULLINATED PEPTIDE; ANTI CCP	1200	Serum	2ml	Rf/Fro	Overnight fasting is preferred.	sample by fri report on sat
61	ANTI-BODY TO Jo-1	1500	Serum	2ml	Rf/Fro		7 Days
62	ANTI-BODY TO Scl-70; SCLERODERMA ANTIBODY	1500	Serum	2ml	Rf/Fro		7 Days
66	ANTI -CETROMERE ANTIBODY, SERUM	1500	Serum	2ml	Rf/Fro		7 days
68	ANTI-DIURETIC HORMONE ADH	4000	WB-EDTA	4ml	Frozen	4 ml plasma from 2 lavender top (EDTA) tube immediately put in ice, centrifuge in a refrigerated centrifuge and separate in transport vials freeze. Ship frozen. DO NOT THAW	8 Days
69	ANTI-ds DNA ANTIBODY (Double Stranded)	750	Serum	2ml	Rf/Fro		Next day
71	ANTI-ENDOMYSIAL ANTIBODY, IgA	1100	Serum	2ml	Rf/Fro	3mL (1.5mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	3 days
72	ANTI-GLIADIN ANTIBODY, IgA	1100	Serum	2ml	Rf/Fro		Next day
73	ANTI-GLIADIN ANTIBODY, IgG	1500	Serum	2ml	Rf/Fro		One week
74	ANTI-GLUTAMIC ACID DECARBOXYLASE ANTIBODIES (GAD-65)	5000	Serum	2ml	Rf/Fro		One week
75	ANTI-HISTONE ANTI BODIES	1500	Serum	2ml	Rf/Fro		One week
76	ANTI-MITOCHONDRIAL ANTIBODY, IMMUNOFLUORESCENCE; AMA, IFA	1100	Serum	2ml	Rf/Fro		7 Days
77	ANTI-NUCLEAR ANTIBODY / FACTOR, EIA; ANA/ANF, SERUM, EIA	550	Serum	2ml	Rf/Fro	Overnight fasting is preferred.	Next Day
78	ANTI-NUCLEAR ANTIBODY / FACTOR, IMMUNOFLUORESCENCE; ANA / ANF, IFA (HEP-2)	1000	Serum	2ml	Rf/Fro	Overnight fasting is preferred.	Sample by Tue/ Fri Report on Wed / Sat
79	ANTI -PARIETAL CELL ANTIBODY (GASTRIC)	1500	Serum	2ml	Rf/Fro	Overnight fasting is preferred.	7 days
80	ANTI-PHOSPHOLIPID ANTIBODY, IgG AND IgM	500 each	Serum	2ml			Next day

# CAREONE DIAGNOSTICS

## SPECIAL TEST PRICE-LIST

81	ANTI-RNP-Sm IgG, Anti ribonucleoprotein smith antigen complex	1500	Serum	2ml	Rf/Fro		One week
82	ANTI-Sm IgG, Anti Smith IgG	1500	Serum	2ml	Rf/Fro		One week
83	ANTI-SMOOTH MUSCLE ANTIBODY, IMMUNOFLUORESCENCE; ASMA, IFA	1100	Serum	2ml	Rf/Fro		Sample by Tue/Fri Reporting Wed/Sat
84	ANTI-SPERM ANTIBODY, SERUM	800	Serum	2ml	Rf/Fro	Overnight fasting is preferred.	Next Day
85	ANTI-STREPTOLYSIN - O; ASLO; ASO TITRE	250	Serum	2ml	Rf/Fro		Same day
86	ANTI-TAENIA SOLIUM, IgG, SERUM; CYSTICERCOSIS ANTIBODY	800	Serum	2ml	Rf/Fro		Next day
87	ANTI -TREPONEMA PALLIDUM, IgG, SERUM, FTA - Abs	2000	Serum	2ml	Rf/Fro		7 days
88	ANTI-THYROGLOBULIN ANTIBODY; ANTI Tg	850	Serum	2ml	Rf/Fro		Same day
89	ANTI-THYROID PEROXIDASE ANTIBODY; ANTI TPO, (ANTI TMA THYROID MICROSOMAL ANTIBODY)	850	Serum	2ml	Rf/Fro		Same day
90	APOLIPOPROTEINS A1 & B Includes Apo A1/Apo B Ratio.	500	Serum	2ml	Rf/Fro	Minimum 12-hour fasting is required.	Next Day
91	Arsenic , blood Serum	2500	WB-EDTA	5ml	Ambient	5 ml (3 ml min.) Royal blue tube available from lab	10 Days
92	Arsenic, 24 hrs urine	2500	Urine	10ml of the Total collection	Rf/Fro	10 ml aliquot of 24 hr urine collected in a special metal free container collected from lab patient should avoid all herbal medications and sea food for 3 days prior to sample collection. Ship refrigerated or frozen	10 Days
93	ASPERGILLUS ANTIBODY IgE SPECIFIC	1500	Serum	2ml	Rf/Fro		7 Days
94	Aspergillus antibody IgG, Serum	1200	Serum	2ml	Rf/Fro		7 Days

# CAREONE DIAGNOSTICS

## SPECIAL TEST PRICE-LIST

95	Aspergillus antibody IgM, Serum	1200	Serum	2ml	Rf/Fro		7 Days
98	BARBITURATES RANDOM URINE	500	Urine	10ml	Rf/Fro	10ml urine sample ( 5 ml min) random fresh voided.	Next day
100	Benzodiazepine Screen, Random urine (Qualitative)	500	Urine	10ml	Refrigerated	Aliquot of random urine sample.	Next Day
101	Benzodiazepine , Urine, Quantative	600	Urine	10ml	Refrigerated	Aliquot of random urine sample.	3 days
102	BETA - 2 - MICROGLOBULIN, SERUM	850	Serum	3ml	Rf/Fro		3 days
108	Blood culture aerobic only (Rapid)	550	Bactec	5ml	Ambient	5 ml blood in aerobic bottle available from lab	prelim 24 hrs, final 5 days
109	Blood culture anerobic (Rapid)	650	Bactec	5ml	Ambient	5 ml blood in anerobic bottle provided by lab	prelim - 24 hrs, final 5 days
111	BLOOD GAS ANALYSIS, ARTERIAL	550	Heparinised Syringe&needle	2ml	Refrigerated	2 mL arterial blood in a preheparin- ised (1000 IU/mL) disposable syringe and needle. Avoid air bubbles. Mix gently 10 times. Seal needle with rubber stopper. Ship refrigerated. DO NOT FREEZE.	Same Day
115	BONE MARROW EXAMINATION	500	EDTA Bone Marrow, Peripheral Blood	2ml/4ml	Ambient	Submit 1 mL (0.5 mL min.) Bone Marrow in 1 Green Top (Sodium Heparin) tube OR minimum 4 smears of Bone marrow aspirate AND 2 peripheral blood smears, 4 ML EDTA BLOOD. Fix smears in 95 % ethanol for 30 min and air dry. Give clinical history on test requesition form	Next day
116	Bone marrow reticuline stain on biopsy	500	Bone Marrow Slide,or Trepchine Biopsy	Block/Ustained Slide	Ambient	wax block of trephine biopsy or , 4 unstained slides , trephine biopsy	72 hrs
119	BRUCELLA ANTIBODY; IgG .	1000	Serum	2ml	Rf/Fro		7 Days

# CAREONE DIAGNOSTICS

## SPECIAL TEST PRICE-LIST

120	BRUCELLA ANTIBODY; IgM	1000	Serum	2ml	Rf/Fro		7 Days
125	CA 15.3; BREAST CANCER MARKER	850	Serum	2ml	Rf/Fro		Same day
127	CA24.2, CANCER MARKER, SERUM	1300	Serum	2ml	Rf/Fro		5 days
128	CA 72.4, GASTRIC CANCER MARKER	1200	Serum	2ml	Rf/Fro		5 days
129	CALCITONIN	1400	Serum	2ml	Rf/Fro	2 mL (0.5 mL min.) serum from 1 Red Top (No Additive) tube. Ship frozen. DO NOT THAW. Overnight fasting is preferred.	4 days
133	CANNABINOIDES, SCREEN RANDOM URINE MARIJUANA, SCREEN, RANDOM URINE	500	Serum	2ml	Refrigerated	Random urine	next day
134	CARBAMAZEPINE LEVELS	500	Serum	2ml	Refrigerated	Provide Last dose, preferd sample before the next dose	Same day
135	CATECHOLAMINES, PLASMA	3000	Plasma-Heparin	8ml	Frozen	8 ML (4 ML min) Plasma from 4 special green top available from lab separate plasma within 10 min of collection. Transfer into plastic tranport container and freeze. Ship frozen. Do not Thaw	5 days
136	CATECHOLAMINES, 24-HOUR URINE Includes Epinephrine, Norepinephrine and Dopamine.	3000	Urine	50ml	Ambient	50mL (10mL min.) aliquot of 24-hour urine. Collect urine with appropriate volume of 50% HCl as per the chart on page 8 on Specimen Collection, to maintain pH between 1-3. Do not use concentrated HCl. Record 24 hour volume on test request form and urine container.	5 days
137	CD4/CD8 /CD3 COUNT	1800	WB-EDTA	5ml	Ambient	5 ml , 2 lavender top (EDTA) tubes, ship at 18-22 degree celsius	4 days
138	CD4/CD8 COUNT	1400	WB-EDTA in 2 Vials	5ml	Ambient	5 ml , 2 lavender top (EDTA) tubes, ship at 18-22 degree celsius	4 days
142	CHLAMYDIA TRACHOMATIS DNA PCR	4000	Any specimen	2ml	Ambient	Any specimen	7 days
143	CHLAMYDIA TRACHOMATIS IgG	900	Serum	2ml	Rf/Fro	Overnight fasting is preferred.	5 days
144	CHLAMYDIA TRACHOMATIS Ig A	900	Serum	2ml	Rf/Fro		5 days



# CAREONE DIAGNOSTICS

## SPECIAL TEST PRICE-LIST

145	CHLAMYDIA TRACHOMATIS Ig M	1500	Serum	2ml	Rf/Fro	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	7 Days
156	CHYLE IN FLUID	500	Fluid	10ml	Rf/Fro	10 ml of fluid in a sterile container and 2 ml serum in 1 SST tube. Ship refrigerated. Do not freeze	3 days
161	CMV AVIDITY TEST	1000	Serum	2ml	Rf/Fro	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	7 Days
164	CMV-DNA PCR-RT	4000	Any specimen	3ml	Ambient	Any Specimen	7 Days
165	COBALT BLOOD	3000	Royal Blue Tube Heparin	5ml	Refrigerated	Whole blood in(Royal blue Top Heparin) tube available from lab	One week
166	COBALT URINE 24 HRS/ RANDOM	3000	Urine	10ml	Refrigerated	10 ml (5 ml min) aliquot of 24 hr urine collected in acid washed metal free container available from lab	One week
167	COCAINE RANDOM URINE	500	Urine			10 ml (5 ml min) aliquot of random urine in a sterile container	Next Day
171	COPPER, 24-HOUR URINE	1000	Urine	10ml	Refrigerated	10 mL (5 mL min.) aliquot of 24-hour urine collected in an acid-washed (metal free) container kit available from lab. Ship refrigerated or frozen.	Next day
172	COPPER, SERUM	1000	Serum	3ml	Refrigerated	3 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	next day
173	CORTISOL, FREE, 24-HOUR URINE	500	Urine	10ml	Refrigerated	10 mL (5 mL min.) aliquot of 24-hour urine. Collect urine with 10 grams of boric acid to maintain sample integrity. Record 24-hour volume on test request form and urine container. Ship refrigerated or frozen.	Next Day
174	CORTISOL, SERUM	350	Urine	2ml	Refrigerated	2 mL (1 mL min.) serum from 1 SST. Diurnal variation present; take sample between 4-6 p.m. for evening sample and 8 - 10 AM for morning sample Specify time and date on sample container and test request form. Ship refrigerated or frozen.	Next day
175	C-PEPTIDE, SERUM	700	Serum	2ml	Frozen	2mL (0.5mL min.) serum from 1 SST. Overnight fasting is required. Ship frozen. DO NOT THAW.	3rd day
182	CRYOGLOBULINS, SCREENING (QUALITATIVE)	600	Serum	3ml	Ambient	3 mL (1.5 mL min.) serum. Collect blood in a PRE-WARMED Red Top tube using a PRE-WARMED syringe and needle at 37° C. Allow to clot and centrifuge preferably at 37° C. Separate serum and transport at	5 days

# CAREONE DIAGNOSTICS

## SPECIAL TEST PRICE-LIST

						room temperature. DO NOT FREEZE	
183	CRYPTOCOCCUS ANTIGEN, CSF	1100	CSF	2ml	Frozen	2 mL (0.5 mL min.) cerebrospinal fluid in a sterile screw capped container. Ship frozen. DO NOT THAW.	same day
184	CRYPTOCOCCUS ANTIGEN, SERUM	1100	Serum	2ml	Frozen	2 mL (0.5 mL min.) serum from 1 Red Top (No Additive) tube. Do not use SST gel barrier tubes. Ship frozen. DO NOT THAW.	Same day
185	CRYPTOSPORIDIUM / ISOSPORA IDENTIFICATION, STOOL Ziehl Neelsen stains.	150	Stool	5gm	Ambient	5 g (1 g min.) of fresh stool in special transport medium available from lab. Ship at 18-22°C.	Next day
186	CRSTALS IN SYNOVIAL FLUID	120	Fluid	2ml	Refrigerated	1 ML fluid in a sterile container, ship refrigerated.	Next day
197	CULTURED AFB (MYCOBACTERIA), ANTIBIOTIC SENSITIVITY, RAPID - 5 Drugs Panel (Streptomycin, Isoniazid, Rifampicin, Ethambutol)	3000				Submit pure growth of AFB on appropriate media in screw capped container in a sealed plastic bag. Ship refrigerated OR, when used with other AFB culture or Mycobacteria identification tests, the pure growth obtained in the Laboratory from specimens already processed	3 weeks
198	CULTURED AFB (MYCOBACTERIA), ANTIBIOTIC SENSITIVITY, RAPID - 10 Drugs Panel (Streptomycin, Isoniazid, Rifampicin, Ethambutol, Kanamycin, Cycloserine, Ofloxacin/ Ciprofloxacin Pas)	6000				Submit pure growth of AFB on appropriate media in screw capped container in a sealed plastic bag. Ship refrigerated OR, when used with other AFB culture or Mycobacteria identification tests, the pure growth obtained in the Laboratory from specimens already processed	3 weeks
199	CYCLOSPORINE, TROUGH LEVEL (CO)	1500	WB-EDTA	3ml	Refrigerated	3 ml (1.5 ml min) whole blood in a lavender top (EDTA) tube ship refrigerated. DO Not Freeze. Recommended specimen just before the next dose	next day
200	CYCLOSPORINE PEAK LEVEL (C2)	1500	WB-EDTA	3ml	Ambient	3 ml (1.5 ml min.) whole blood recommended sample 2 hr after the dose at steady state.	next Day
201	CYFRA 21-1 LUNG CANCER MARKER	2000	Serum	2ml	Ambient		7 Days
207	CSF PROTEIN ELECTROPHORESIS OLIGOCCONAL BAND	3500	CSF	2ml		CSF in a sterile container.	7 Days
210	DHEA ; DEHYDROEPIANDRO- STERONE ( DHEA), SERUM	2500	Serum	2ml		2 mL (0.5 mL min.) serum from 1 SST. Overnight fasting is preferred. Specify age and sex on test request form. Ship refrigerated or frozen.	One week

# CAREONE DIAGNOSTICS

## SPECIAL TEST PRICE-LIST

212	DIAZEPAM (BENZODIAZEPINE) QUALITATIVE, URINE, RANDOM	500	Urine	10ml		10 ml ( 5 ml, min) randomurine aliquot.	Next Day
213	DIGOXIN (Lanoxin)	600	Serum	2ml		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Do not use SSTgel barrier tubes. Separate serum within 2 hours of collection. Ship refrigerated or frozen. Trough: just before the next dose is due at steady state. Peak: 6-8 hours after dose	3 days
214	DI- HYDROTESTOSTERONE (DHT)	1200	Serum	2ml		2 mL (0.5 mL min.) serum from 1 SST. Overnight fasting is preferred. Specify age and sex on test request form. Ship refrigerated or frozen.	One week
215	DILANTIN SODIUM/ (PHENYTOIN) / EPTOIN, SERUM	550	Serum	2ml		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Do not use SSTgel barrier tubes. Separate serum within 2 hours of collection. Ship refrigerated or frozen. Trough: just before the next dose is due at steady state. Peak: 6-8 hours after dose	Same day
218	ECHOVIRUS IgG	2000	Serum	3ml	Rf/Fro		2 weeks
219	ECHOVIRUS IgM	2000	Serum	3ml	Rf/Fro		2 weeks
220	ECHINOCOCCUS IgG, SERUM	800	Serum	3ml	Rf/Fro		Next Day
223	EPO, ERYTHROPOIETIN	1300	Serum	3ml	Rf/Fro		4 days
224	EPSTEIN-BARR VIRUS ANTIBODY IgG	1200	Serum	3ml	Rf/Fro		10 Days
225	EPSTEIN-BARR VIRUS ANTIBODY IgM	1200	Serum	3ml	Rf/Fro		10 Days
226	EPTOIN (PHENYTOIN), SERUM	550	Serum	3ml	Rf/Fro		Same day
228	ESTRIOL, TOTAL, SERUM; E3, TOTAL, SERUM	600	Serum	3ml	Rf/Fro	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen	3 days
229	ESTRIOL, UNCONJUGATED, SERUM; uE3,	1000	Serum	3ml	Rf/Fro	2 mL (1 mL min.) serum from 1 SST. Ship	3 days

# CAREONE DIAGNOSTICS

## SPECIAL TEST PRICE-LIST

	SERUM					refrigerated or frozen.	
230	ESTROGEN RECEPTOR; ER .	1000	Parafin Block			Submit tumour tissue in 10% Formal saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology/Cytopath	Block/Sections/Smears : 5 days Tissues : 7 days
231	ETHANOL (QUANTITATIVE)	500	Serum	2ml	Refrigerated	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen	Same day
232	FACTOR XIII, CLOT SOLUBILITY	250	Plasma-Citratd	4ml	Frozen	4 mL (3 mL min.) Platelet Poor Plasma (PPP) from 2 Blue Top (Sodium Citrate) tubes prepared within 1 hour of collection. Aliquot PPP in 2 transport vials and FREEZE. Ship frozen. DO NOT THAW. See page 6 for instructions on PPP preparation. Overnight fastin	next day
233	FACTOR IX, FUNCTIONAL	2500	Plasma-Citratd	4ml	Frozen	4 mL (3 mL min.) Platelet Poor Plasma (PPP) from 2 Blue Top (Sodium Citrate) tubes prepared within 1 hour of collection. Aliquot PPP in 2 transport vials and FREEZE. Ship frozen. DO NOT THAW. See page 6 for instructions on PPP preparation. Overnight fastin	1 week
234	FACTOR VIII FUNCTIONAL	2500	Plasma-Citratd	4ml	Frozen	4 mL (3 mL min.) Platelet Poor Plasma (PPP) from 2 Blue Top (Sodium Citrate) tubes prepared within 1 hour of collection. Aliquot PPP in 2 transport vials and FREEZE. Ship frozen. DO NOT THAW. See page 6 for instructions on PPP preparation. Overnight fastin	1 week
240	FLOW CYTOMETRY, INDIVIDUAL MARKER - CD 3	1300	WB-EDTA,WB-HEPARIN	3ml	Ambient	3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL ( 2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2mL (1mL min.) heparinized Bone Marrow. Ship immediately at 18-22°C. DO NOT FREEZE. Specify time, date and clinical detail	3 days
241	FLOW CYTOMETRY, INDIVIDUAL MARKER - CD 4	1400	WB-EDTA,WB-HEPARIN	3ml	Ambient	3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL ( 2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2mL (1mL min.) heparinized Bone Marrow. Ship	3 days

# CAREONE DIAGNOSTICS

## SPECIAL TEST PRICE-LIST

						immediately at 18-22°C. DO NOT FREEZE. Specify time, date and clinical detail	
242	FLOW CYTOMETRY, INDIVIDUAL MARKER - CD 8	1300	WB-EDTA,WB-HEPARIN	3ml	Ambient	3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2mL (1mL min.) heparinized Bone Marrow. Ship immediately at 18-22°C. DO NOT FREEZE. Specify time, date and clinical detail	3 days
243	FLOW CYTOMETRY, CD3/CD4/CD8	1800	WB-EDTA,WB-HEPARIN	3ml	Ambient	3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2mL (1mL min.) heparinized Bone Marrow. Ship immediately at 18-22°C. DO NOT FREEZE. Specify time, date and clinical detail	3 days
244	FLOW CYTOMETRY, INDIVIDUAL MARKER - HLA B27	1500	WB-EDTA,WB-HEPARIN	3ml	Ambient	3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2mL (1mL min.) heparinized Bone Marrow. Ship immediately at 18-22°C. DO NOT FREEZE. Specify time, date and clinical detail	3 days
251	FOLATE , RBC; FOLIC ACID RBC	2000	WB-EDTA,Serum	2ml/1ml	Refrigerated	2 mL (1 mL min.) whole blood from 1 Lavender Top (EDTA) tube AND 2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Wrap tubes in aluminium foil to protect from light. Ship refrigerated. DO NOT FREEZE. Record PCV on test request form.	Next day
252	FOLATE, SERUM; FOLIC ACID, SERUM	750	Serum	2ml	Refrigerated	2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Wrap tube in aluminium foil to protect from light. Ship refrigerated or frozen.	Same day
253	FREE BETA HCG, SERUM	900	Serum	3ml	Rf/Fro	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	3 days
254	FREE PSA, SERUM	750	Serum	3ml	Rf/Fro	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	3 days
257	FREE TESTOSTERONE	900	Serum	3ml	Rf/Fro	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	7 days
261	FTA- ANTIBODIES (IgG ANTIBODIES TO TREPONEMA)	2000	Serum	3ml	Rf/Fro	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	10 Days

# CAREONE DIAGNOSTICS

## SPECIAL TEST PRICE-LIST

263	GALL STONE ANALYSIS	1000	Stone		Ambient	Submit stone (s) in a screw capped plastic container. Ship at room temperature.	7 days
264	GASTRIN	1200	Serum	3ml	Rf/Fro	2 mL (1 mL min.) serum from 1 SST. Ship frozen. DO NOT THAW. Overnight fasting is required.	5 days
275	GLUCOSE-6-PHOSPHATE DEHYDROGENASE, QUANTITATIVE; G-6-PD, QUANTITATIVE	600	WB-EDTA	3ml	Refrigerated	3 mL (1.5 mL min.) whole blood in 1 LavenderTop (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Next day
277	GOMORI'S METHENAMINE SILVER STAIN	500	FNAC	Block/Ustained Slide	Ambient	Submit tissue biopsy in 10% formal saline OR formalin fixed paraffin embedded tissue block OR FNAC specimen smeared onto 4 glass slides and air dried. Ship at room temperature. Brief clinical history in Histopathology Requisition Form is mandatory.	4 days
279	GROWTH HORMONE; GH; hGH; HUMAN GROWTH HORMONE	500	Serum	3ml	Rf/Fro	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Next day
286	HAEMOSIDERIN, URINE	200	Urine	10ml	Ambient	10 mL (5 mL min.) aliquot of random urine. First morning specimen is preferred. Ship refrigerated. DO NOT FREEZE. OR Centrifuge urine at 1200g for 15 min. Smear deposit on four glass slides, air dry, fix in 95% ethanol for 30 min & ship at room temperature	Next day
287	HAM TEST; PAROXYSMAL NOCTURNAL HEMOGLOBINURIA, CONFIRMATORY TEST; PNH CONFIRMATORY TEST;	1000	Serum & WB-EDTA	3ml	Refrigerated	3 mL (1.5 mL min.) serum from 1 SST AND 3 mL (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	3 days
288	HAPTOGLOBULIN, SERUM	1000	Serum	3ml	Rf/Fro	2ML(1ML MIN) serum from 1 SST, Ship refrigerated or frozen	3 days
294	HELICOBACTER PYLORI IgG ANTIBODY	1000	Serum	3ml	Rf/Fro		6 days
295	HELICOBACTER PYLORI IgA ANTIBODY	1000	Serum	3ml	Rf/Fro		7 days
296	HEPATITIS A, TOTAL	700	Serum	3ml	Rf/Fro		3 days
297	HEPATITIS A IgM ANTIBODY; ANTI- HAV, IgM	700	Serum	3ml	Rf/Fro		Same day
298	HEPATITIS B CORE, IgM ANTIBODY; HB CORE IgM; HBc IgM	700	Serum	3ml	Rf/Fro		Same day
299	HEPATITIS B CORE, TOTAL ANTIBODY; HB CORE TOTAL; HBc TOTAL	700	Serum	3ml	Rf/Fro		Same day
300	HEPATITIS B SURFACE ANTIBODY; ANTI HBs; HBsAb	500	Serum	3ml	Rf/Fro		Next day

# CAREONE DIAGNOSTICS

## SPECIAL TEST PRICE-LIST

301	HEPATITIS B SURFACE ANTIGEN; HBsAg; AUSTRALIA ANTIGEN Screening only.	200	Serum	3ml	Rf/Fro		Same day
302	HEPATITIS Be ANTIBODY; ANTI - HBe; HBe Ab	550	Serum	3ml	Rf/Fro		Next day
303	HEPATITIS Be ANTIGEN; HBe Ag	550	Serum	3ml	Rf/Fro		Same day
305	HEPATITIS C ANTIBODY, IgG; ANTI - HCV, IgG	1000	Serum	3ml	Rf/Fro		5 days
306	HEPATITIS C ANTIBODY, IgM; ANTI HCV, IgM	1000	Serum	3ml	Rf/Fro		5 days
307	HEPATITIS E ANTIBODY, IgG; ANTI-HEV, IgG; HEV IgG	1500	Serum	3ml	Rf/Fro		6 days
308	HEPATITIS E ANTIBODY, IgM; ANTI-HEV, IgM; HEV IgM	1000	Serum	3ml	Rf/Fro		Next day
315	HSV - DNA Detction for (1+2)	3000	Serum	3ml	Ambient	Serum	7 days
325	HIV 1 & 2 ANTIBODIES, WESTERN BLOT	2200	Serum	2ml	Refrigerated	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	4 Days
327	HIV- p24 ANTIGEN	1100	Serum	2ml	Refrigerated	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	3 days
328	HIV - Proviral DNA	2500	WB-EDTA	5ml	Refrigerated	EDTA Blood	5 days
329	HLA - B27, DISEASE ASSOCIATION , FLOW CYTOMETRY	1500	WB-EDTA/WB-HEPARIN	4ml	Ambient	4 mL (2.5 mL min.) whole blood in 1 Lavender Top (EDTA/Green Top (Sodium Heparin) tube OR 10 mL (5 mL min.) whole blood in 1 Yellow Top (ACD) tube. Ship at 18-22°C. DO NOT FREEZE.	3 days
330	HOMOCYSTEINE, SERUM, QUANTITATIVE	700	Serum	2ml	Refrigerated	2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen.	Same day
331	HOMOGENITISIC ACID (ALKEPTANURIA), URINE	600	Urine	10ml	Refrigerated	10 ML of fress voided urine in a sterile container, wwrap in black paper or aluminium foil to avoid exposure to light., ship refrigerated	2nd day
332	IGF -1 , INSULIN LIKE GROWTH (SOMATOMEDIN C)	3000	Serum	2ml	Frozen	2 mL (1 mL min) serum from 1 SST. Separate serum in a transport vial and freez, Ship frozen. DO NOT THAW	3 days
333	IGF BINDING PROTEIN- 3	3500	Serum	2ml	Frozen	2 mL (1 mL min) serum from 1 SST. Separate serum in a transport vial and freez, Ship frozen. DO NOT THAW	3 days
334	IMMUNOFIXATION ELECTROPHORESIS, 24-HOUR URINE	6500	Urine	50ml	Refrigerated	50 mL (30 mL min.) aliquot of 24-hour urine without preservative. Record 24 hour volume on test request form and urine container. Ship refrigerated or frozen.	7 days

# CAREONE DIAGNOSTICS

## SPECIAL TEST PRICE-LIST

335	IMMUNOFIXATION ELECTROPHORESIS, SERUM	5500	Serum	2ml	Refrigerated	2 mL (1 mL min.) serum from 1 SST. Separate serum within 1 hour of collection. Ship refrigerated or frozen. Grossly hemolysed specimens are not accepted.	7 days
336	IMMUNOGLOBULIN IgA, SERUM	400	Serum	2ml	Refrigerated	2 mL (1 mL min.) serum from 1 SST. Overnight fasting is preferred. Ship refrigerated or frozen.	3 days
338	IMMUNOGLOBULIN IgG, SERUM	400	Serum	3ml	Rf/Fro	2 mL (1 mL min.) serum from 1 SST. Overnight fasting is preferred. Ship refrigerated or frozen.	3 days
339	IMMUNOGLOBULIN IgG CSF	800	CSF	2ml	Refrigerated	2 ML CSF in a sterile container, ship refrigerated or frozen.	3 days
340	IMMUNOGLOBULIN IgM, SERUM	400	Serum	3ml	Rf/Fro	2 mL (1 mL min.) serum from 1 SST. Overnight fasting is preferred. Ship refrigerated or frozen.	3 days
342	INTRINSIC FACTOR ANTIBODY	2000	Serum	3ml	Rf/Fro	2 mL (1 mL min.) serum from 1 SST. Ship frozen. DO NOT THAW	6 days
343	INFECTIOUS MONONUCLEOSIS ANTIBODY; MONOSPOT TEST Qualitative test for heterophile antibodies to Infectious Mononucleosis	700	Serum	3ml	Rf/Fro	2 mL (1 mL min.) serum from 1 SST. Ship frozen. DO NOT THAW	Next day
347	INSULIN ANTIBODIES	2000	Serum	3ml	Rf/Fro	2 mL (1 mL min) serum from 1 SST. Ship refrigerated or frozen.	7 days
349	ISLET CELL ANTIBODY BY IMMUNOFLOURESCENCE	2500	Serum	3ml	Rf/Fro	2 mL (1 mL min) serum from 1 SST. Ship refrigerated or frozen.	7 days
351	KIDNEY STONE ANALYSIS ; RENAL/URETERIC CALCULUS ANALYSIS	800	Stone		Refrigerated	Submit stone (s) in a clean screw capped container. Ship at room temperature or refrigerated.	7 Days
354	LA ANTIBODEIS ( SSB- SOLUTE SUBSTANCE B)	1500	Serum	2ml	Refrigerated	2 mL (1 mL min) serum from 1 SST. Ship refrigerated or frozen.	5 days
359	LEAD, 24 HRS URINE	1500	Urine	10ml	Refrigerated	10 ML ( 5 ML MIN) urine aliquot of 24 hrs urine collected in apsecial metal free container provided by lab. Ship refrigerated or frozen.	1 week
360	LEAD, BLOOD	1500	WB-ROYAL BLUE Vial	3ml	Refrigerated	2ML( 1ML MIN) Whole blood from aspecial metal free royal blue top heparine tube available from lab.Ship refrigerated. DO NOT FREEZE	1 Week
361	LEISHMANIA ANTIBODY IgG; KALA AZAR ANTIBODY	1200	Serum	3ml	Rf/Fro	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	One week
362	LEPTOSPIRA ANTIBODY, IgG	900	Serum	3ml	Rf/Fro	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	3 days



# CAREONE DIAGNOSTICS

## SPECIAL TEST PRICE-LIST

363	LEPTOSPIRA ANTIBODY, IgM	900	Serum	3ml	Rf/Fro	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	3 days
371	LIPID PROFILE, EXTENDED 1 Includes Lipoprotein (a) Apolipoproteins A1, B & Ratio *Cholesterol *Triglycerides *HDL/VLDL/LDL *Cholesterol/HDL Ratio *Chylomicrons,	1500	Serum	3ml	Refrigerated	3 mL (2 mL min.) serum from 1 SST. Separate serum within 1 hour of collection. Ship refrigerated or frozen. Minimum 12 hours overnight fasting is mandatory. Grossly hemolysed samples are not acceptable.	3 days
372	LIPID PROFILE DIRECT includes direct estimation of HDL, LDL, Cholesterol, triglyceride	650	Serum	3ml	Refrigerated	3 mL (2 mL min.) serum from 1 SST. Separate serum within 1 hour of collection. Ship refrigerated or frozen. Minimum 12 hours overnight fasting is mandatory. Grossly hemolysed samples are not acceptable.	Next day
373	LIPOPROTEIN (a); Lp(a)	900	Serum	3ml	Rf/Fro	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Same day
376	LIVER KIDNEY MICROSOMAL ANTI BODY, IMMUNOFLOUORESCENCE; LKM ANTIBODY	1100	Serum	3ml	Rf/Fro		1 Week
377	LUPUS ANTICOAGULANT PANEL; LAC PANEL DIL. APTT, BY DRRVT, BY KCT	1600	Plasma-Citratd	6ml	Frozen	6 mL (4.5 mL min.) Platelet Poor Plasma (PPP) from 3 Blue Top (Sodium Citrate) tubes prepared within 1 hour of collection. Aliquot PPP in 4 transport vials and FREEZE. Ship frozen. DO NOT THAW. See page 6 for instructions on PPP preparation. Overnight fas	Next Day
378	MAGNESIUM, RANDOM URINE	500	Urine	10ml	Refrigerated	10 mL (5 mL min.) aliquot of 24-hour urine in a plastic screw capped container. Do not use preservatives. Record 24-hour volume on test request form and urine container. Ship refrigerated or frozen.	3 days
379	MAGNESIUM, 24-HOUR URINE	500	Urine	10ml	Refrigerated	10 mL (5 mL min.) aliquot of 24-hour urine in a plastic screw capped container. Do not use preservatives. Record 24-hour volume on test request form and urine container. Ship refrigerated or frozen.	3 days
380	MAGNESIUM, SERUM	250	Serum	3ml	Rf/Fro		Same day
385	MARIJUANA (CANNABIS )	500	Urine	10ml	Ambient		next day
386	METANEPHRINES, 24-HOUR URINE *Metanephriene * Normetanephriene	2500				50 mL (10 mL min.) aliquot of 24-hour urine. Collect urine with appropriate volume of 50% HCl as per the chart on page 8 on Specimen Collection, to maintain pH between 1-3. Do not use concentrated HCl. Record 24 hour volume on test request	5 days

# CAREONE DIAGNOSTICS

## SPECIAL TEST PRICE-LIST

						form and urine	
387	MEASLES (RUBEOLA) - IgG ANTIBODIES, Serum	1200	Serum	3ml	Rf/Fro		One week
388	MEASLES (RUBEOLA) - IgM ANTIBODIES, Serum	1200	Serum	3ml	Rf/Fro		One week
389	MEASLES (RUBEOLA) - IgM ANTIBODIES, CSF	1200	Csf	2ml	Refrigerated		One week
390	MEASLES (RUBEOLA) - IgG ANTIBODIES, CSF	1200	Csf	2ml	Refrigerated		One week
394	MICROSOMAL ANTIBODIES ( ANTI TPO)	850	Serum	3ml	Rf/Fro		Next day
395	MITOCHONDRIAL ANTIBODIES	1000	Serum	3ml	Rf/Fro		7 Days
396	MONOSPOT TEST , INFECTIOUS MONONUCLEOSIS	700	Serum	3ml	Rf/Fro		Next Day
397	MUCICARMINE STAIN, HISTOPATHOLOGY	500				Submit tissue biopsy in 10% formal saline OR formalin fixed paraffin embedded tissue block OR FNAC specimen smeared onto 4 glass slides and air dried. Ship at room temperature. Brief clinical history in Histopathology Requisition Form is mandatory.	4 days
398	MUMPS VIRUS IgG ANTIBODY	1200	Serum	3ml	Rf/Fro		One week
399	MUMPS VIRUS IgM ANTIBODY	1200	Serum	3ml	Rf/Fro		One week
404	OPIATES, (MORPHINE)	500	Urine	10ml	Ambient		Next day
405	OSMOLALITY, RANDOM URINE	600	Urine	5ml	Refrigerated	5 mL (3 mL min.) aliquot of random urine. DO NOT USE PRESERVATIVE. Ship refrigerated only. DO NOT FREEZE.	Next day
406	OSMOLALITY, SERUM	600	Serum	2ml	Refrigerated	2 mL (1 mL min.) serum from 1 SST. Separate serum within 1 hour of collection. Ship refrigerated or frozen. Grossly hemolysed samples not accepted.	Next day
407	OSTEOCALCIN - SERUM	2000	Plasma-Heparin	3ml	Frozen	3 ml plasma from a green top heparine tube, separate plasma immediately and ship frozen.	3 days
408	OXALATE, 24-HOUR URINE	2000	Urine	10ml	Refrigerated	10 mL (5 mL min.) aliquot of 24 hour urine. Collect urine with 25 mL of 50% HCl to maintain pH below 3. For children adjust the volume of acid according to the expected urinary output. Keep urine	10 Days

# CAREONE DIAGNOSTICS

## SPECIAL TEST PRICE-LIST

						refrigerated during collection. Mix well before aliquoting.	
413	PAP SMEAR WITH PROCEDURE	300				Patient to report to lab OR one of its Satellite Labs for the test.	Next Day
414	PARATHYROID HORMONE ( INTACT PTH)	1000	Serum	3ml	Rf/Fro	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	3 days
415	PARIETAL CELL ANTIBODY, TOTAL	1200	Serum	3ml	Rf/Fro	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	7 Days
418	PHENCYCLIDINE (PCP), RANDOM URINE	500	Urine	10ml	Refrigerated	10 ml aliquot of random urine sample, ship refrigerated or frozen	same day
419	PHENOBARBITONE (Luminal, Gardenal )	550	Serum	2ml	Refrigerated	2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Do not use SST gel barrier tubes. Ship refrigerated or frozen. Recommended specimen: just before the next dose is due at steady state.	Same day
420	PHENYTOIN (Dilantin, Epilam plain, Epileptin, Eptoin, Epsolin)	550	Serum	2ml	Refrigerated	2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Do not use SST gel barrier tubes. Ship refrigerated or frozen. Recommended specimen: just before the next dose is due at steady state.	Same day
425	PLATELET ANTIBODY	3000	Serum	2ml	Frozen	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	10 Days
426	PORPHOBILINOGEN, QUALITATIVE, URINE, 2nd MORNING SAMPLE; PBG - Screen	1000	Urine	10ml	Refrigerated	Collect 2nd morning 10mL (5 mL min.) urine sample. Aliquot into a dark coloured vial/container wrapped completely with black paper or aluminium foil, containing 0.1 mL Glacial Acetic Acid. Ship refrigerated or frozen.	2 days
427	PORPHYRINS, 24 HRS URINE - Quantitative	4000	Urine	10ml	Refrigerated	10 ml aliquot of 24 well mixed urine without preservative collected and transported in a dark coloured container or container wrapped in aluminium foil or black paper	7 Days
436	PROTEIN C- FUNCTIONAL ASSAY	3500	Plasma-Citrated	4ml	Frozen	4 ml ( 2 ml min) platelets poor plasma from 2 blue top citrate tube. Separate plasma within 1 hr of collection in transport vials and freeze. Ship frozen. DO NOT THAW.	8 days
437	PROTEIN S FUNCTIONAL ASSAY	3500	Plasma-Citrated	4ml	Frozen	4 ml ( 2 ml min) platelets poor plasma from 2 blue top citrate tube. Separate plasma within 1 hr of collection in transport vials and freeze. Ship frozen. DO NOT THAW.	8 days
438	PROTEIN ELECTROPHORESIS, 24-HOUR URINE Does not include IFE.	3500	Urine	20ml	Refrigerated	20 mL (10 mL min.) aliquot of well mixed 24 - hour urine collection without preservative.	One week

# CAREONE DIAGNOSTICS

## SPECIAL TEST PRICE-LIST

						Record total volume on test request form and urine container. Ship refrigerated only. DO NOT FREEZE.	
439	PROTEIN ELECTROPHORESIS, SERUM Includes total protein, electrophoretic fractions with densitometric tracing. Does not include IFE.	700	Serum	2ml	Refrigerated	2 mL (1 mL min.) serum from 1 SST. Separate serum within 1 hour of collection. Ship refrigerated only. DO NOT FREEZE. Grossly hemolysed samples not accepted.	4 days
445	PSA, FREE ; PROSTATE - SPECIFIC ANTIGEN, FREE Prostatic Cancer Marker. Along with total PSA +RATIO OF TWO	1000	Serum	2ml	Frozen	2 mL (1 mL min.) serum from 1 SST. Ship frozen. DO NOT THAW Do not draw sample within 7 days of Digital Rectal Examination (DRE) or Rectal Prostatic Ultrasonography. .	3 days
447	PTH, INTACT; PARATHYROID HORMONE, INTACT	1000	Serum	3ml	Frozen	3 mL (1.5 mL min.) serum from 1 SST between 8-10 a.m Allow blood to clot. Separate serum and freeze immediately within 30 minutes. Ship frozen. DO NOT THAW. Hemolysed samples are not acceptable.	Next day
448	RBC FRAGILITY TEST; OSMOTIC FRAGILITY TEST FOR RBC-S	500	WB-Heparin ,Peripheral Smears	3ml	Ambient	3 mL (1.5 mL min.) whole blood from 1 Green Top (Sodium Heparin) tube AND 2 peripheral blood smears. Ship refrigerated. DO NOT FREEZE.	2 days
451	RETICULIN STAIN, HISTOPATHOLOGY	500	Tissue Biopsy, Formalin Fixed Tissue		Ambient	Submit tissue biopsy in 10% formal saline OR formalin fixed paraffin embedded tissue block OR FNAC specimen smeared onto 4 glass slides and air dried. Ship at room temperature. Brief clinical history in Histopathology Requisition Form is mandatory.	4 days
457	RUBELLA AVIDITY TEST	1000	Serum	3ml	Rf/Fro		4 days
462	SEX HORMONE BINDING GLOBULIN; SHBG	1300	Serum	3ml	Rf/Fro	Specify age and sex on test request form.	3 Days
466	SPERM ANTIBODY (ASAB) TOTAL	1000	Serum/Semen	2ml/1ml	Refrigerated		Next day
474	SSA ( SOLUBLE SUBSTANCE - A) ANTI BODY	1500	Serum	3ml	Rf/Fro		7 Days
475	SS B ( SOLUBLE SUBSTANCE - B) ANTI BODY	1500	Serum	3ml	Rf/Fro		7 Days
482	SUCROSE LYSIS TEST; PAROXYSMAL NOCTURNAL HEMOGLOBINURIA SCREENING TEST; PNH SCREENING TEST;	500	WB-EDTA	3ml	Refrigerated	DO NOT FREEZE.	2 days

# CAREONE DIAGNOSTICS

## SPECIAL TEST PRICE-LIST

488	TESTOSTERONE, FREE	1000	Serum	3ml	Rf/Fro	Specify age and sex on test request form.	3 days
489	TESTOSTERONE, TOTAL	400	Serum	2ml	Rf/Fro	Specify age and sex on test request form.	2nd day
491	THROMBIN TIME; TT	900	WB-Citrate	3ml	Frozen	For Local samples: 3mL (2mL min.) whole blood in 1 Blue Top (Sodium Citrate) tube with proper dilution. Ship at 18-22°C, to reach lab within 3 hours. For Outstation samples: 2mL (1.5mL min.) Platelet Poor Plasma (PPP) from 2 Blue Top (Sodium Citrate) tube separate PPP in transport vials and freeze. Ship frozen. DO NOT THAW	Next day
492	THYROGLOBULIN; Tg Cancer Marker	1500	Serum	2ml	Frozen	DO NOT THAW. At least 6 weeks should elapse after thyroidectomy or needle biopsy, or I131 therapy before Specimen collection.	Next day
493	THYROGLOBULIN ANTIBODIES	850	Serum	3ml	Rf/Fro		Same day
495	TISSUE TRANSGLUTAMINASE ANTIBODY, IgA; tTG Antibody, IgA	1000	Serum	3ml	Rf/Fro		next day
502	TOXOPLASMA AVIDITY	1000	Serum	3ml	Rf/Fro		5 days
503	TRANSFERRIN	700	Serum	3ml	Rf/Fro	Overnight fasting is preferred.	Same day
504	TREPONEMA PALLIDUM HEMAGGLUTINATION; TPHA	350	Serum	3ml	Rf/Fro		4 days
509	TROPONIN - I	1100	Serum	3ml	Frozen	DO NOT THAW.	Same Day
510	TROPONIN - T	1100	Serum	3ml	Frozen	DO NOT THAW.	Same Day
517	UROBILINOGEN, QUALITATIVE, 24-HOUR URINE	250	Urine	10ml	Rf/Fro	10 mL (5 mL min.) aliquot of 24-hour urine without preservative, in a dark coloured container. Ship refrigerated or frozen. Keep Specimen refrigerated during collection and protect from light.	Same day
518	VALPROIC ACID; VALPROATE (Epilex, Valparin, Encorate, Depakene)	500	Serum	2ml	Ambient	Do not use SST gel barrier tubes. Recommended Specimen: just before the next dose is due at steady state.	Same day
519	VARICELLA ZOSTER VIRUS IgG ANTIBODY; VZV IgG ANTIBODY (Chicken Pox; Herpes zoster; Shingles)	1500	Serum	3ml	Rf/Fro		One week
520	VARICELLA ZOSTER VIRUS IgM ANTIBODY; VZV IgM ANTIBODY (Chicken Pox; Herpes zoster; Shingles)	1500	Serum	3ml	Rf/Fro		One week
523	VITAMIN A	2000	Serum	2ml	Rf/Fro	PROTECT FROM LIGHT BY WRAPPING IN BLACK PAPER OR ALUMINIUM FOIL.	One week

# CAREONE DIAGNOSTICS

## SPECIAL TEST PRICE-LIST

524	VITAMIN B12, SERUM; CYANOCOBALAMIN, SERUM	750	Serum	3ml	Rf/Fro	Wrap in aluminium foil to protect from light. Overnight fasting is required.	One week
525	VITAMIN B1, SERUM	2500	Plasma- EDTA	2ml	Rf/Fro	Separate plasma immediately and freeze. DO NOT THAW. PROTECT FROM LIGHT BY WRAPPING IN ALUMINIUM FOIL.	One week
526	VITAMINE B2 (RIBOFLAVIN)	2500	Plasma- EDTA	2ml	Rf/Fro	Separate plasma immediately and freeze. DO NOT THAW. PROTECT FROM LIGHT BY WRAPPING IN ALUMINIUM FOIL.	One week
527	VITAMINE B 6( PYRIDOXINE)	2500	Plasma- EDTA	2ml	Rf/Fro	Separate plasma immediately and freeze. DO NOT THAW. PROTECT FROM LIGHT BY WRAPPING IN ALUMINIUM FOIL.	One week
528	VITAMINE C	2500	Plasma- EDTA	2ml	Rf/Fro	Separate plasma immediately and freeze. DO NOT THAW. PROTECT FROM LIGHT BY WRAPPING IN ALUMINIUM FOIL.	One week
529	VITAMIN D, 1, 25-DIHYDROXY	3500	Serum	3ml	Rf/Fro	Wrap in aluminium foil to protect from light. Overnight fasting is required.	10 Days
530	VITAMIN D3 (25 HYDROXY)	2500	Serum	3ml	Rf/Fro	Wrap in aluminium foil to protect from light. Overnight fasting is required.	One week
531	VITAMIN E (TOCOPHEROL)	2500	Serum	3ml	Rf/Fro	Wrap in aluminium foil to protect from light. Overnight fasting is required.	8 days
532	VMA, 24-HOUR URINE; VANILMANDELIC ACID	2500	Urine, Aliquot of 24-hour urine	50ml	Rf/Fro	Collect urine with appropriate volume of 50% HCl as per the chart on page 8 on Specimen Collection, to maintain pH between 1-3. Do not use concentrated HCl. Record 24 hour volume on test request form and urine	3 days
533	VMA, RANDOM URINE; VANILMANDELIC ACID	2500	Urine	50ml	Rf/Fro	50 mL (10 mL min.) aliquot of random urine. Add 50% HCl drop by drop to bring the pH between 1-3. Do not use concentrated HCl. Ship refrigerated or frozen. Also see list of drugs to be avoided on page 8 on Specimen Collection.	3 days
535	XYLOSE- D 5 HOUR URINE	2000	Urine(Sodium Fluoride)	10ml	Refrigerated	Overnight fasting is required. Patient must void and discard the urine prior to Xylose -D administration. For adults administer orally 25 gm of D - Xylose dissolved in 250 ml of water. Additional 250 ml of water can be given. For children 0.5 gm/kg body weight of D- Xylose can be administered with a maximum 25 gms. Collect all urine voided over next 5 hrs. Send 5 ml aliquot in a grey top tube. Ship refrigerated. Do not	Sample by Tue/Report by Sat

# CAREONE DIAGNOSTICS

## SPECIAL TEST PRICE-LIST

						freeze. Record 5 hr urine volume on test request form.	
536	XYLOSE- D PLASMA	2000	Plasma-Sodium Fluoride	3ml	Refrigerated	Overnight fasting is required. . For adults administer orally 25 gm of D - Xylose dissolved in 250 ml of water. Additional 250 ml of water can be given. For children 0.5 gm/kg body weight of D- Xylose can be administered with a maximum 25 gms. Submit 2 ml plasma from 1 grey top tube (sodium fluoride) taken 2 hrs afetr the orall administration of D- Xylose.. separate plasma within 1 hr of collection. Ship refrigerated DO NOT FREEZE.	Sample by Tue/Report by Sat
537	ZINC SERUM	1100	Serum	2ml	Rf/Fro	Special royal blue top tube availble from lab.	7 days

**PLEASE FEEL FREE TO CALL OUR HELPDESK FOR ANY CLARIFICATION / INQUIRY.**

**HELPDESK:0612-2533690**

**9334832890,9431055234**